



LAB	CODE	TEST NAME	CPT4 CODE
ACL	ZIKAPCR	Zika Virus by PCR, Blood Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 6 Weeks	87662
ACL	ZIKA AB	Zika Virus IgM Antibody Capture (MAC) by ELISA Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year (avoid repeated freeze/thaw cycles)	86794
ACL	ZPP	Zinc Protoporphyrin Container :One lavender (EDTA) 3.0 mL (also acceptable: one royal blue (trace elements, EDTA) 6.0 mL) Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: Unacceptable	84202
ACL	ZINC 8	Zinc Transporter 8 Antibody Container :One gold gel 3.5 mL (also acceptable: one plain red 6.0 mL) Stability :Ambient: 5 Days / Refrigerated: 1 Week / Frozen: 3 Months	86341
ACL	ZIN	Zinc, Blood Container :One royal blue (trace elements, EDTA) 6.0 mL (also acceptable: one royal blue (trace elements, no additive) 6.0 mL or one dark blue (sodium heparin) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 2 Weeks / Frozen: 2 Weeks	84630
ACL	ZNRBC	Zinc, RBC Container :One royal blue (trace elements, EDTA) 6.0 mL Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: Unacceptable	84630
ACL	ZON	Zonisamide Level Container :One plain red 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 1 Week / Refrigerated: 1 Week / Frozen: 6 Weeks	80203