

Bronchial washing/Bronchoalveolar lavage (BAL)

Each specimen container needs to have a *minimum* of two (2) patient identifiers to be accepted.

- Label clean leak-proof specimen container with patient's full name (first and last name), date of birth or medical record number (MRN). Include specimen source on container.
- Place at least 2 mL of fresh specimen in container and tighten cap.
- Place in refrigerator until transport to lab.
- An accompanying requisition with patient's full name (first and last), date
 of birth, medical record number (MRN), specimen source, ordering
 physician's name, date of collection, diagnosis, and pertinent patient
 history.
- If Oil Red O for lipid-laden macrophages is desired, please indicate in the comments field in the EPIC Non-GYN order or on the accompanying requisition form. Specimen must be <u>fresh</u> without fixative added.
- If iron stain for hemosiderin-laden macrophages is desired please indicate in the comments field in the EPIC Non-GYN order or on the accompanying requisition form. Specimen must be <u>fresh</u> without fixative added.
- If Grocott's methenamine silver stain (GMS) on cell block is desired, please indicate in the comments field in the EPIC Non-GYN order or on the accompanying requisition form. Specimen must be <u>fresh</u> without fixative added.