Legacy Laboratory Services

Legacy LabAlert

December 2023

An Important Update from Legacy Laboratory Services

Specimen Transport Bag

Effective December 20th, 2023, specimen bags with Legacy logos have been replaced with Labcorp standard bags. Please see below for pictures of the new specimen bags and use their corresponding PeopleSoft numbers (PS#) going forward.

If you have Legacy-branded bags, continue to use them.

Please continue to use the Legacy-labeled lockbox for all Legacy-directed samples.

Table 1: Specimen Bag Information

Bag picture	Description	Ordering information	Quantity per order
INACAGE SECTION ONLY INACAGE SECTION ONLY	Bag 6X10 <u>ROOM</u> <u>TEMP</u> barcoded LCA	PS #: 129862	50/pack
SPECIMEN STRICT REFRIGERATE LabCorp	Bag 6X10 REFRIGERATE 3-wall tearzone LCA	PS #: 110475	100/pack



 $Save\ a\ copy\ to\ \underline{https://legacyhealth.sharepoint.com/sites/\underline{LaboratoryServicesHub/Lab}\ Alerts/Forms/\underline{AllItems.aspx}}$

Page 1 of 2 Revision: 1.9.2023

Legacy Laboratory Services

SPECIMEN ONLY FROZEN	Bag 6X10 <u>FROZEN</u> 3- wall tearzone LCA	PS #: 19930	100/pack

Order Form Links:

Legacy Clinics and Labs

Independent Clinics

For additional information, please contact your account representative, client services or consult the test table on our website.

Legacy Laboratory Client Services: 503-413-1234, 877-270-5566, www.legacyhealth.org\labservices



Page 2 of 2 Revision: 9.24.20