

OnLine Test Menu



CODE	TEST NAME	CPT4 CODE
LAB550	OBSTETRIC PANEL Container:One 6 ml EDTA (pink), one 4 ml EDTA (lavender), one 5 ml Sserum separator tube (gold) Performing Laboratory:CHI Health Laboratory	
LAB5710	OBSTETRIC PANEL WITH HIV Container:One 6 ml EDTA (pink), one 4 ml EDTA (lavender), one 5 ml Sserum separator tube (gold) Performing Laboratory:CHI Health Laboratory	
LAB696	OCCULT BLOOD GASTRIC / DUODENUM Performing Laboratory:CHI Health Laboratory	82273
LAB694	OCCULT BLOOD X 1, STOOL Performing Laboratory:CHI Health Laboratory	82270
LAB3447	OLIGOCLONAL BAND PROFILE Performing Laboratory:LABCORP	
LAB3147	OPIATES CONFIRMATION URINE Performing Laboratory:LABCORP SPOKANE	
LAB418	ORGANIC ACIDS URINE Container:Sterile container Temperature:CRITICAL Frozen. Separate specimens must be submitted when multiple tests are ordered. Performing Laboratory:ARUP	83918
LAB107	OSMOLALITY, SERUM Container:One 5 ml serum separator tube (Gold) or one 3 ml plasma separator tube (light green) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	83930
LAB420	OSMOLALITY, URINE Container:Sterile container Performing Laboratory:CHI Health Laboratory	83935
LAB1134	OSMOTIC FRAGILITY, RBC Container:lavender(EDTA) Temperature:Refrigerated Performing Laboratory:ARUP	85555
LAB955	OVA AND PARASITES EXAMINATION Container:Stool in SAF (fill container to appropriate line) Temperature:Room temp Performing Laboratory:CHI Health Laboratory	87177 88313
LAB421	OXALATE, URINE, 24 HOUR Container:24 hour urine container Temperature:CRITICAL FROZEN Performing Laboratory:CHI Health Laboratory	83945
LAB484	OXCARBAZEPINE LEVEL Performing Laboratory:LABCORP BURLINGTON	