

CODE	TEST NAME	CPT4 CODE
LAB548	SACCHAROMYCES CEREVISIAE ANTIBODIES, IGG AND IGA Performing Laboratory:LABCORP BURLINGTON	
LAB34	SALICYLATE LEVEL Container:One 5 mL serum separator tube (gold) or one 3 mL plasma separator tube (light green) (Min: 1.0 mL whole blood) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	80329
LAB1235	SCHISTOSOMA ANTIBODY, IGG Performing Laboratory:CAMBRIDGE BIOMEDICAL INC	
LAB322	SEDIMENTATION RATE, AUTOMATED Container:EDTA (lavender) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	85652
LAB579	SELENIUM SERUM Container:Royal blue (No Additive), Royal blue (K2EDTA), or Royal Blue (Na2EDTA) Due to nationwide shortages of some phlebotomy tubes specified for trace element testing, ARUP Laboratories is currently accepting red top (serum) and Lavender (EDTA) tubes. Temperature:Room temperature. Also acceptable: refrigerated or frozen. Performing Laboratory:ARUP	84255
LAB216	SEMEN ANALYSIS Performing Laboratory:Physicians Laboratory	
LAB891	SEMEN ANALYSIS, POST-VASECTOMY Performing Laboratory:Physicians Laboratory	
LAB4225	SERINE PROTEASE3, IGG Performing Laboratory:LABCORP BURLINGTON	
LAB121	SEROTONIN SERUM Performing Laboratory:LABCORP BURLINGTON	
LAB4016	SEX HORMONE BINDING GLOBULIN Performing Laboratory:LABCORP DENVER	
LAB1318	SHIGA-LIKE TOXIN ANTIGEN, EIA Container:Para-Pak Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	87427
LAB339	SICKLE CELL SCREEN Container:EDTA (lavender) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	85660
LAB344	SJOGRENS SYNDROME-A EXTRACTABLE NUCLEAR ANTIBODY Container:5 ml serum separator tube (Gold) or 3 ml plasma separator tube (Light green) Performing Laboratory:CHI Health Laboratory	86235
LAB345	SJOGRENS SYNDROME-B EXTRACTABLE NUCLEAR ANTIBODY Container:5 ml serum separator tube (Gold) or 3 ml plasma separator tube (Light green) Performing Laboratory:CHI Health Laboratory	86235
LAB2364	SMITH (ENA) AB, IGG Container:5 ml serum separator tube (Gold) or 3 ml plasma separator tube (Light green) Performing Laboratory:CHI Health Laboratory	86235