

LAB	CODE	TEST NAME
Evergreen	GABA-LC	Gabapentin Level
Evergreen	GGT	Gamma Glutamyl Transferase Container:Gold/Tiger (serum) Alternate: Mint green (lithium heparin); Red (serum)
Evergreen	FGQ1B	Ganglioside GQ1b Ab (IgG) Container:Red Top Preferred/SST Temperature:Refrigerate
Evergreen	GAST-LC	Gastrin Level
Evergreen	GENTPK	Gentamicin Level, Peak Container:Red (serum) Alternate: Mint green (lithium heparin); Gold/Tiger (serum)and other gel-type tubes, however, they may artificially, randomly lower results if they are not promptly centrifuged and separated.
Evergreen	GENTR	Gentamicin Level, Random Container:Red (serum) Alternate: Mint green (lithium heparin); Gold/Tiger (serum)
Evergreen	GENTTR	Gentamicin Level, Trough Container:Red (serum) Alternate: Lithium heparin plasma (green top tube), SST and other gel-type tubes, however, they may artifactually, randomly lower results if they are not promptly centrifuged and separated.
Evergreen	GLUM	Gestational Glucose 1Hr Scrn Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Gray top tube (NaF-KOx); Mint green (lithium heparin)
Evergreen	GLAG	Giardia Antigen Container:Cary-Blair or sterile leak proof container. Temperature:Store and transport as appropriate for particular specimen type.
Evergreen	GLBMAB-LC	Glomerular Basement Membrane Antibody
Evergreen	LAB82	GLUCOSE RANDOM Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Gray top tube (NaF-KOx); Mint green (lithium heparin) Temperature:Store and transport refrigerated.
Evergreen	GLU2PP	Glucose Tolerance 2 Hour (non-pregnant) Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Gray top tube (NaF-KOx); Mint green (lithium heparin) Temperature:Store and transport refrigerated.
Evergreen	3 Hour Glucose (non-pregnant)	Glucose Tolerance, 3 HR (non-pregnant) Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Mint Green (lithium heparin)
Evergreen	LAB169	Glucose Tolerance, Pregnant (2HR) Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Gray top tube (NaF-KOx); Mint green (lithium heparin)
Evergreen	GEST3	Glucose Tolerance, Pregnant (3HR) Container:Gold/Tiger (serum) or Gray top (NaF-KOx) Alternate: Gray top tube (NaF-KOx); Mint green (lithium heparin)