

| LAB       | CODE      | TEST NAME                                                                                                                                                                                       |
|-----------|-----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Evergreen | HPBTAD-LC | H Pylori Breath Test                                                                                                                                                                            |
| Evergreen | HAPT-LC   | Haptoglobin                                                                                                                                                                                     |
| Evergreen | HCGUR     | HCG Qualitative Urine<br>Container: Sterile leakproof plastic urine container<br>Temperature: Store and transport refrigerated or freeze (-20) if transport will exceed 2 days.                 |
| Evergreen | HCGQL     | HCG Qualitative, Serum<br>Container: Gold/Tiger (serum), Red (serum)<br>Temperature: Store and transport refrigerated or freeze (-20) if transport will exceed 2 days.                          |
| Evergreen | HCGQT     | HCG Quantitative, Serum<br>Container: Gold/Tiger (serum), Red (serum), Mint Green (lithium heparin)<br>Temperature: Store and transport refrigerated or freeze if transport will exceed 3 days. |
| Evergreen | HVYMTL-LC | Heavy Metals Panel, WB                                                                                                                                                                          |
| Evergreen | HSSUB-LC  | Helper Suppressor                                                                                                                                                                               |
| Evergreen | HCT       | Hematocrit<br>Container: Lavender (EDTA, 4mL)                                                                                                                                                   |
| Evergreen | HGB       | Hemoglobin<br>Container: Lavender (EDTA, 4mL)                                                                                                                                                   |
| Evergreen | HGBF-LC   | Hemoglobin F Quant                                                                                                                                                                              |
| Evergreen | HGBTha-LC | Hemoglobinopathy Fractionation Panel                                                                                                                                                            |
| Evergreen | HEMO      | Hemogram Only<br>Container: Lavender (EDTA, 4mL)                                                                                                                                                |
| Evergreen | LAB20     | HEPATIC FUNCTION PANEL<br>Container: Gold/Tiger (serum)    Alternate: Mint green (lithium heparin); Red (serum)                                                                                 |
| Evergreen | HAAB-LC   | Hepatitis A Antibody                                                                                                                                                                            |
| Evergreen | HAABM-LC  | Hepatitis A IgM                                                                                                                                                                                 |