

LAB	CODE	TEST NAME
Evergreen	NTX-LC	N Teloepptides, Urine
Evergreen	LAB13999	NASH FibroSure(R) Plus
Evergreen	NMS1	Necrotizing Myopathy Evaluation, S Container:Red Top Preferred/SST Temperature:Refrigerate
Evergreen	NEOCBC	Neonatal CBC Container:Lavender Top Microtainer (EDTA) 7ml EDTA tubes are NOT preferred; can cause a delay in testing.
Evergreen	LAB13391	Neurofilament Light Chain
Evergreen	PKU	NEWBORN SCREENING (WASHINGTON) / PKU Container:PKU card Temperature:Room temperature
Evergreen	Nicotine-LC	Nicotine Metabolite, Urine
Evergreen	NMO-LC	NMO IgG Auto Ab
Evergreen	NMRLIP-LC	NMR Lipoprofile Test
Evergreen	BNP	NT-proBNP Container:Preferred: Mint green (lithium heparin) Alternate: Lavender EDTA Plasma or Serum
LABCORP	183624	N. gonorrhoeae NAA, Confirm Container:Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) or SurePath(TM) liquid cytology vials Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)
LABCORP	205914	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.
LABCORP	222000	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.
LABCORP	831560	N-Acetyl-L-Aspartic Acid,Urine Container:Plastic urine container Temperature:FREEZE
LABCORP	701195	NALOXONE,UR,MS CONFIRM Container:Plastic urine container without preservative Temperature:Room Temperature