

## Online Test Menu Laboratory Catalog



| LAB | 0005   | TEST NAME  | ODT4 00DE |
|-----|--------|--|-----------|
| LAB | CODE   | TEST NAME  | CPT4 CODE |
| GMC | HPYC   | H Pylori Culture w/reflex<br>Place Of Service:LabCorp  |           |
| GMC | HFLUBG | H. influenzae b Antibody, IgG<br>Turn around time:2 - 4 days<br>Place Of Service:LABCORP BURLINGTON    |           |
| GMC | HPYBR  | H. Pylori Breath Test<br>Place Of Service:LabCorp  |           |
| GMC | HANTA  | Hantavirus Antibodies<br>Turn around time:6 - 7 days<br>Place Of Service:QUEST DIAG. INFECT. DISEASE   |           |
| GMC | HAPT   | Haptoglobin Turn around time:1 day Place Of Service:LABCORP SEATTLE                                    |           |
| GMC | HCGQ   | HCG, QUANTITATIVE Turn around time:NEXT DAY Place Of Service:GRITMAN MEDICAL CENTER                    | 84702     |
| GMC | HCGTM  | HCG, Quantitative Tumor Marker Turn around time:within 1 day Place Of Service:DYNACARE NORTHWEST INC   |           |
| GMC | HCVFIS | HCV Fibrosure Turn around time:3 - 5 days Place Of Service:LABCORP BURLINGTON                          |           |
| GMC | HCVGTY | HCV Genotype by PCR&Probe(Reflex) Turn around time:3 - 5 days Place Of Service:LABCORP BURLINGTON      |           |
| GMC | HCVPGT | HCV RNA Quant by PCR Reflex Geno<br>Turn around time:3 - 5 days<br>Place Of Service:LABCORP BURLINGTON |           |
| GMC | HDLD   | HDL, DIRECT Turn around time:NEXT DAY Place Of Service:GRITMAN MEDICAL CENTER                          | 83718     |
| GMC | HMGHGB | HEALTHY MEASURES GHGB Turn around time:NEXT DAY Place Of Service:GRITMAN MEDICAL CENTER                | 83036     |
| GMC | HMPNL  | HEALTHY MEASURES PANEL Turn around time:NEXT DAY Place Of Service:GRITMAN MEDICAL CENTER               | 80061     |
| GMC | UHMR4  | Heavy Metals 4, Urine Random<br>Turn around time:3 - 6 days<br>Place Of Service:LABCORP SPOKANE        |           |
| GMC | HVYMTL | Heavy Metals Panel, Whole Bld<br>Turn around time:2 - 5 days<br>Place Of Service:LABCORP SPOKANE       |           |