

LAB	CODE	TEST NAME	CPT4 CODE
KSB	EOSNA	Nasal Eosinophil Container:Nasal Pharangeal Swab	89190
KSB	NEOPT	Neopterin, Serum - LC	
KSB	NEO	Neoscreen Panel Container:PKU paper Temperature:Room Temp	84030
KSB	UNICO	Nicotine & Cotinine, Urine - LC	
KSB	NICO	Nicotine and Metabolite, Quant - LC	
KSB	NORO	Norovirus, RT-PCR - LC	
KSB	NTELU	N-Telopeptide, Urine - LC	
LABCORP	183624	N. gonorrhoeae NAA, Confirm Container:Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) or SurePath(TM) liquid cytology vials Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87591
LABCORP	247315	N.Meningitidis A/C/W-135/Y IgG Container:Gel-barrier tube or red-top tube Temperature:Refrigerated	86741 x4
LABCORP	222000	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84300 84133 82436
LABCORP	205914	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84300 84133 82436
LABCORP	810960	Nabumetone as Metabolite Container:Red-top tube, lavender-top (EDTA) tube fluoride) tube Temperature:Refrigerated Stability: Room Temperature: 30 days Refrigerated: 30 days Frozen: 17 months	80299
LABCORP	831560	N-Acetyl-L-Aspartic Acid,Urine Container:Plastic urine container Temperature:FREEZE	82542
LABCORP	805241	NailStat-12 Container:Contact lab for Hair collection kit. Temperature:Maintain specimen at room temperature.	80307
LABCORP	805131	NailStat-14 Container:Contact lab for Hair collection kit. Temperature:Maintain specimen at room temperature.	80307