

LAB	CODE	TEST NAME
MSHS	70212	LACOSAMIDE
MSHS	6934	LACTALBUMIN, ALPHA IGE Container:Gold Top Tube Temperature:Refrigerate
MSHS	1146	LACTATE Container:•Sodium heparin tube •Sodium fluoride/potassium oxalate tube Temperature:On Ice or 2-8 / Analyze fresh sample.
MSHS	3221	LACTATE CSF Container:Sterile Container Temperature:Frozen
MSHS	74035	Lactate CSF
MSHS	72226	LACTIC ACID, CSF
MSHS	70170	LACTOFERRIN, FECAL, QN
MSHS	8593	LACTOFERRIN, QUAL, STOOL Container:Leak proof sterile plastic screw top container without additive Temperature:Refrigerate
MSHS	5308	LAMB/MUTTON Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.
MSHS	5307	LAMB/MUTTON -F88 Container:Gold Top Tube Temperature:Refrigerate
MSHS	70137	LAMOTRIGINE, SERUM
MSHS	5437	LATEX ENHANCED Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.
MSHS	1444	LD FLUID Container:Red Top Tube Temperature:Refrigerate
MSHS	70516	LD ISOENZYMES
MSHS	1130	LD(LDH) Container:Gold Top Serum Separator Tube, Plain Red Top Tube (No Additive) Temperature:Refrigerate