

LAB	CODE	TEST NAME
MSHS	5446	OAK (WHITE) Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.
MSHS	5445	OAK (WHITE) -T7 Container:Gold Top Tube Temperature:Refrigerate
MSHS	5332	OAT Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.
MSHS	5331	OAT -F7 Container:Gold Top Tube Temperature:Refrigerate
MSHS	70619	OCCULT BLOOD, FECAL, IA
MSHS	2040	OCCULT BLOOD, STOOL Container:Hemocult Slide Temperature:Cards should be sent to lab ASAP
MSHS	70596	OLANZAPINE (ZYPREXA)
MSHS	70538	OLIGOCLONAL BANDING
MSHS	60345	ONCOLOGY ARRAY Container:Bone Marrow (CYTO GRN w/ MEDIA+HEPARIN)/ Peripheral Blood (Green Top Sodium Heparin Vacutainer) / Fresh Tumor (CYTO GRN w/ MEDIA+HEPARIN) / Spleen (Container with Phosphate-Buffered-Saline (PBS) free of Calcium and Magnesium) Temperature:Room Temperature or 37C
MSHS	5335	ONION Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.
MSHS	5334	ONION -F48 Container:Gold Top Tube Temperature:Refrigerate
MSHS	1204	OPIATE, URINE Container:Urine Container Temperature:Refrigerate
MSHS	52024	OPIATES CONFIRMATION, URINE
MSHS	70597	OPIATES, UNCONJUGATED EXPANDED
MSHS	5338	ORANGE Container:Gold-top (Serum tubes (with or without gel barrier) Temperature:Analyze fresh as possible.