

LAB	CODE	TEST NAME
MSHS	2519	V.W.F. ANTIGEN Container:Blue Top Tube Temperature:Refrigerate
MSHS	2518	V.W.F. RISTO:COF Container:Blue Top Tube Temperature:Specimen must be transported to the laboratory within 4 hrs of draw at room temperature
MSHS	70621	VAGINITIS/VAGINOSIS, DNA PROBE
MSHS	1270	VALPROIC ACID Container:Serum tubes (with or without gel barrier) Temperature:2-8 / Analyze fresh sample.
MSHS	74616	VALPROIC ACID (TOTAL+FREE)
MSHS	1199	VANCOMYCIN TROUGH Container:Serum Tubes (with or without gel barrier or clot activator) Temperature:2-8 / Analyze fresh sample.
MSHS	70372	VANILLYLMANDELIC ACID, 24-HR U
MSHS	70020	VARICELLA ZOSTER IGG
MSHS	70101	VARICELLA-ZOSTER AB, IGM
MSHS	73414	VARICELLA-ZOSTER VIRUS (VZV) AB (TOTAL, IGM), ACIF/IFA, CSF
MSHS	6891	VARICELLA-ZOSTER VIRUS DNA Container:Sterile Container Temperature:Refrigerate
MSHS	5039	VARICELLA-ZOSTER VIRUS DNA Container:Sterile Container Temperature:Frozen
MSHS	424	VDRL (CSF ONLY) Container:Leak proof sterile plastic screw top container without additive Temperature:Refrigerate
MSHS	70116	VEGF, PLASMA
MSHS	72210	VENLAFAXINE (EFFEXOR)