



Online Test Menu

Laboratory Catalog

LAB	CODE	TEST NAME	CPT4 CODE
Mission Health	P2Y12	P2Y12 Inhibitor Container:Contact lab for special collection packet.	85576
Mission Health	PNELTLC	Pancreatic Elastase, Fecal LC Place Of Service:LabCorp	
Mission Health	PAPSMEAR	Pap Smear Screen Container:1 each ThinPrep (obtain collection container from Lab Pathology)	88142
Mission Health	PAPSMEARA	Pap Smear, Anal Container:1 each ThinPrep (obtain collection container from Lab Pathology)	88175
Mission Health	PARASLC	Paraneoplastic Ab LC Place Of Service:LabCorp	
Mission Health	PARAP	Parasite Panel	87328 87329
Mission Health	PTHW	Parathyroid Hormone on Vitros 7600 XT, 2LAB.CHM.5138.00 Container:4 mL SST or 4 mL Green Lithium Heparin Temperature:RT - ≤1 day REF - ≤10 DAYS Frozen - ≤6MTHS	83970 82310
Mission Health	PTHIO	Parathyroid Hormone, Intraoperative Container:3 mL Lavender EDTA Tube	83970
Mission Health	PRVLC	Parvovirus B19, Human, IgG/IgM LC Place Of Service:LabCorp	
Mission Health	PATHTISS	Pathology Surgical Pathology Specimens Container:Surgical Pathology Specimens 1) Collect tissue specimen. 2) Place tissue sample into a container of formalin. Tissue to formalin ratio should be 1:10. 3) Close lid completely to prevent specimen and fluid from leaking. 4) Label container with patient name, second identifier and specimen source. 5) Place container into biohazard bag and seal. 6) Complete requisition and place into specimen bag pouch.	
Mission Health	PBGRLC	PBG Deaminase, RBC LC Place Of Service:LabCorp	82657
Mission Health	BPH	pH, Body Fluid Container:1 mL of body fluid in Lithium Heparin Tube	83986
Mission Health	FPH	pH, Fecal Container:Collect 1 mL of FRESH stool (no preservatives). Transport to the laboratory in a Sterile Specimen Container.	83986
Mission Health	PHSTLC	pH, Stool LC	
Mission Health	UPH	pH, Urine	81003