

LAB	CODE	TEST NAME	CPT4 CODE
PAMC	LAB581	ZINC Place of Service:LABCORP SPOKANE	
PAMC	LAB1085	ZINC PROTOPORPHYRIN Place of Service:LABCORP SPOKANE	
PAMC	LAB504	ZONISAMIDE LEVEL Place of Service:LABCORP BURLINGTON	
LABCORP	790660	Zaleplon, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.	80299
LABCORP	790710	Zaleplon, Urine Container:Plastic urine container without preservative	80368
LABCORP	489435	ZAP-70 Panel (ASR) Container:Green-top (heparin) tube OR lavender-top (EDTA) tube	88187 88184 88185 x4
LABCORP	163084	Zika Virus Antibody, IgM Container:Gel-barrier tube OR Red-Top tube OR Serum Transfer tube	86794
LABCORP	010180	Zinc Protoporphyrin (ZPP), Workplace Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube.	84202
LABCORP	070111	Zinc Random Urine	84630 82570
LABCORP	070029	Zinc, RBC Container:Royal blue-top (EDTA) tube	84630
LABCORP	003434	Zinc, Urine Container:Plastic urine container, no preservative	84630 82570
LABCORP	070032	Zinc, Whole Blood Container:Royal blue-top (EDTA) tube	84630
LABCORP	790384	Ziprasidone, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.	80299
LABCORP	503995	ZNT8 Antibodies Container:Red-top tube, gel-barrier tube or green-top(heparin) tube	86341
LABCORP	808135	Zolpidem, Serum or Plasma Container:Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended.	80299