

PMG Laboratory Spokane

Test Directory

LAB	CODE	TEST NAME	CPT4 CODE
PMGLAB	LAB239	C BACT URINE Container: Sterile leakproof plastic urine container and then transfer to a urine boric acid tube Temperature: Ambient (room temperature) or refrigerated (unpreserved specimens must be refrigerated) Place Of Service: SHMC	
PMGLAB	LAB7080	C. TRACHOMATIS AND N. GONORRHOEAE, NAAT Container: BD Molecular Swab Collection Kit/BD Molecular Urine Transport Kit Temperature: Room temperature or refrigerate Place Of Service: PMG-PCS	87591 87491 87661
PMGLAB	LAB1376	C. TRACHOMATIS AND N. GONORRHOEAE, NAAT Container: BD Molecular Swab Collection Kit/BD Molecular Urine Transport Kit Temperature: Room temperature or refrigerate Place Of Service: PMG-PCS	
PMGLAB	LAB202	C. TRACHOMATIS AND N. GONORRHOEAE, NAAT (APTIMA) Place Of Service: LABCORP SEATTLE	
PMGLAB	LAB23862	C. TRACHOMATIS AND N. GONORRHOEAE, NAAT Pharyngeal Container: BD MAX UVE Specimen Collection Kit Temperature: Refrigerate Place Of Service: PMG-PCS	87591 87491
PMGLAB	LAB24660	C. TRACHOMATIS AND N. GONORRHOEAE, NAAT RECTAL Container: BD MAX UVE Specimen Collection Kit Temperature: Refrigerate Place Of Service: PMP-PCS	87591 87491
PMGLAB	LAB16058	C1 ESTERASE INHIBITOR TOTAL Place Of Service: LABCORP BURLINGTON	
PMGLAB	LAB769	C1 ESTERASE INHIBITOR, FUNCTIONAL Place Of Service: LabCorp	
PMGLAB	LAB155	CA 125, QUANT Place Of Service: LABCORP SEATTLE	
PMGLAB	LAB776	CA 15-3, QUANT Place Of Service: LabCorp	
PMGLAB	LAB777	CA 19-9, QUANT Place Of Service: LABCORP SEATTLE	
PMGLAB	LAB853	CA 27-29, QUANT Place Of Service: LABCORP SEATTLE	
PMGLAB	LAB519	CALCITONIN Place Of Service: LabCorp	
PMGLAB	LAB53	CALCIUM Container: Red-top tube; SST (gel barrier tube); PST green-top (lithium heparin) Temperature: Refrigerate Place Of Service: PMG-PMP	82310
PMGLAB	LAB54	CALCIUM, IONIZED Container: SST tube, FILLED COMPLETELY (or for infants, 1mL whole blood in a balanced blood gas syringe, on ice) Temperature: Refrigerated. Do not ship on dry ice. Place Of Service: SHMC	