

# PMG Laboratory Spokane

## Test Directory

| LAB    | CODE     | TEST NAME  | CPT4 CODE   |
|--------|----------|--|-------------|
| PMGLAB | LAB25514 | EDOXABAN LEVEL<br>Place Of Service:LabCorp   |             |
| PMGLAB | LAB25368 | EGFR ONCOLOGY, FISH<br>Place Of Service:LabCorp  |             |
| PMGLAB | LAB16    | ELECTROLYTE PANEL<br>Container:Red top tube; SST gel barrier tube; PST gel barrier tube (green lithium heparin)<br>Temperature:Refrigerate<br>Place Of Service:PMG-PMP | 80051       |
| PMGLAB | LAB774   | ENDOMYSIAL AB, IGA<br>Place Of Service:LabCorp   |             |
| PMGLAB | LAB1333  | ENTEROVIRUS, RNA, NAAT<br>Place Of Service:LabCorp   |             |
| PMGLAB | LAB298   | EOSINOPHIL COUNT<br>Container:EDTA lavender top tube<br>Temperature:Refrigerate<br>Place Of Service:PMG-PMP  | 85048       |
| PMGLAB | LAB1812  | EPSTEIN-BARR VIRUS AB, IGG AND IGM<br>Container:SST, red top tube<br>Temperature:Refrigerate<br>Place Of Service:PMG-PCS   | 86664 86665 |
| PMGLAB | LAB1730  | EPSTEIN-BARR VIRUS CAPSID AB, IGG<br>Place Of Service:LABCORP SEATTLE  |             |
| PMGLAB | LAB1731  | EPSTEIN-BARR VIRUS CAPSID AB, IGM<br>Place Of Service:LABCORP SEATTLE  |             |
| PMGLAB | LAB654   | EPSTEIN-BARR VIRUS EARLY ANTIGEN AB, IGG<br>Place Of Service:LABCORP SEATTLE   |             |
| PMGLAB | LAB793   | EPSTEIN-BARR VIRUS NUCLEAR ANTIGEN AB, IGG<br>Place Of Service:LABCORP SEATTLE   |             |
| PMGLAB | LAB863   | EPSTEIN-BARR VIRUS PANEL<br>Place Of Service:LABCORP SEATTLE   |             |
| PMGLAB | LAB2372  | EPSTEIN-BARR VIRUS, NAAT, QUAL<br>Place Of Service:LABCORP BURLINGTON  |             |
| PMGLAB | LAB1373  | EPSTEIN-BARR VIRUS, NAAT, QUANT<br>Place Of Service:LABCORP BURLINGTON   |             |
| PMGLAB | LAB24978 | EPSTEIN-BARR VIRUS, PLASMA, QUANT PCR<br>Place Of Service:LabCorp  |             |