

PMG Laboratory Spokane

Test Directory

LAB	CODE	TEST NAME	CPT4 CODE
PMGLAB	LAB20435	IMMUNOFIXATION, SERUM Place Of Service:LABCORP	
PMGLAB	LAB13407	IMMUNOFIXATION, URINE Place Of Service:LabCorp	
PMGLAB	LAB73	IMMUNOGLOBULIN A Container:PST gel barrier tube (green lithium heparin) Temperature:Refrigerate Place Of Service:PMG-PMP	82784
PMGLAB	LAB74	IMMUNOGLOBULIN E Place Of Service:LABCORP BURLINGTON	
PMGLAB	LAB71	IMMUNOGLOBULIN G Container:PST gel barrier tube (green lithium heparin) Temperature:Refrigerate Place Of Service:PMG-PMP	82784
PMGLAB	LAB12988	IMMUNOGLOBULIN G INDEX, CSF AND SERUM Place Of Service:LABCORP BURLINGTON	
PMGLAB	LAB1001	IMMUNOGLOBULIN G SUBCLASSES Place Of Service:LABCORP SPOKANE	
PMGLAB	LAB72	IMMUNOGLOBULIN M Container:PST gel barrier tube (green lithium heparin) Temperature:Refrigerate Place Of Service:PMG-PMP	82784
PMGLAB	LAB2128	IMMUNOGLOBULIN, FREE LIGHT CHAIN Place Of Service:LabCorp	
PMGLAB	LAB166	IMMUNOGLOBULIN, PANEL, IGG AND IGM AND IGA Container:PST gel barrier tube (green lithium heparin) Temperature:Refrigerate Place Of Service:PMG-PMP	82784 x3
PMGLAB	LAB272	INFLUENZA A AND B AG, IA Container:Dry Nasopharyngeal or Nasal Swab provided only in Immunocard STAT! Flu A&B test kit OR NP/Nares swab in BD Universal Transport Medium Temperature:Refrigerate. Place Of Service:PMG-PMP	
PMGLAB	LAB12992	INFLUENZA A AND B RNA, NAAT Container:Viral transport media Universal Temperature:Refrigerate Place Of Service:PMG-PMP	87502
PMGLAB	LAB25371	INFLUENZA A,B, COVID-19, NAAT Container:BD or Copan Universal Transport Medium (UTM-RT), Remel M4RT, 0.9% physiological saline Temperature:Refrigerate Place Of Service:PMG-PMP	87636
PMGLAB	LAB13217	INHIBIN B Place Of Service:LabCorp	
PMGLAB	LAB649	INSULIN AB Place Of Service:LabCorp	