

PMG Laboratory Spokane

Test Directory

LAB	CODE	TEST NAME	CPT4 CODE
PMGLAB	LAB24783	PAIN MANAGEMENT PANEL 4 REFLEX Place Of Service:LabCorp	
PMGLAB	LAB979	PANCREATIC ELASTASE, FECAL Place Of Service:LABCORP BURLINGTON	
PMGLAB	LAB12010	PAP, LIQUID BASED Container:ThinPrep® or SurePath® vial Temperature:Ambient (room temperature) Place Of Service:SHMC	
PMGLAB	LAB108	PARATHYROID HORMONE, INTACT Container:PST gel barrier tube (green lithium heparin) Temperature:Refrigerate Place Of Service:PMG-PMP	83970
PMGLAB	LAB813	PARATHYROID HORMONE, INTACT AND CALCIUM Container:PST gel barrier tube (green lithium heparin); SST gel barrier tube; Red-top tube Temperature:Refrigerate Place Of Service:PMG LAB-PMP	83970 82310
PMGLAB	LAB800	PARVOVIRUS B19 AB, IGG AND IGM Place Of Service:LABCORP BURLINGTON	
PMGLAB	LAB24688	PE U+IFE U+KAPPA/LAMBDA QN U Place Of Service:LabCorp	
PMGLAB	LAB109	PH, STOOL Place Of Service:LABCORP SEATTLE	
PMGLAB	LAB12489	PH, URINE Container:Leakproof plastic urine container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	
PMGLAB	LAB25462	Phosphatidylethanol (LCorp) Place Of Service:LABCORP	
PMGLAB	LAB1180	PHOSPHATIDYLSERINE AB Place Of Service:LabCorp	
PMGLAB	LAB113	PHOSPHORUS Container:PST gel barrier tube (green lithium heparin); SST gel barrier tube; 1 Full PST or SST microtainer Place Of Service:PMG-PMP	84100
PMGLAB	LAB1177	PLATELET AB, DIRECT Place Of Service:LABCORP	
PMGLAB	LAB301	PLATELET COUNT Container:Lavender top (EDTA) tube Temperature:Room temperature: 24 hours; Refrigerated: 2 days; Frozen: Unacceptable Place Of Service:PMG-PMP	85049
PMGLAB	LAB25369	PML/RARA Fish panel (LCorp) Place Of Service:LABCORP	