

PMG Laboratory Spokane

Test Directory

LAB	CODE	TEST NAME	CPT4 CODE
PMGLAB	LAB2414	T3 REVERSE Place Of Service:LABCORP BURLINGTON	
PMGLAB	LAB137	T3, FREE Container:PST gel barrier tube (green lithium heparin); Red top tube; SST gel barrier tube Temperature:Refrigerate Place Of Service:PMG-PMP	84481
PMGLAB	LAB2168	T3, TOTAL Container:PST gel barrier tube; green top (lithium heparin); SST gel barrier tube; Red top tube Temperature:Refrigerate Place Of Service:PMG-PMP	84480
PMGLAB	LAB23480	T3, TOTAL AND FREE Container:PST gel barrier tube; green top (lithium heparin); SST gel barrier tube; Red top tube Temperature:Refrigerate Place Of Service:PMG-PMP	
PMGLAB	LAB135	T3, UPTAKE Place Of Service:LABCORP SEATTLE	
PMGLAB	LAB127	T4, FREE Container:PST gel barrier tube (green lithium heparin); Red top tube; SST gel barrier tube Temperature:Refrigerate Place Of Service:PMG LAB-KRML	84439
PMGLAB	LAB126	T4, TOTAL Place Of Service:LABCORP SEATTLE	
PMGLAB	LAB876	TACROLIMUS FK506 Place Of Service:LabCorp	
PMGLAB	LAB12323	TACROLIMUS TROUGH LC-MS/MS Place Of Service:LABCORP SEATTLE	
PMGLAB	LAB2280	TBNK LYMPHOCYTE SUBSET Place Of Service:LabCorp	
PMGLAB	LAB343	T-CELL CD4, CD8 COUNT, WITH RATIO Place Of Service:LABCORP	
PMGLAB	LAB23803	TESTOSTERONE TOTAL AND FREE, LC/MS/MS Place Of Service:LABCORP	
PMGLAB	LAB7123	TESTOSTERONE, FREE MS DIALYSIS Place Of Service:ESOTERIX ENDOCRINOLOGY	
PMGLAB	LAB13415	TESTOSTERONE, FREE, DIRECT Place Of Service:LABCORP	
PMGLAB	LAB124	TESTOSTERONE, TOTAL Container:PST gel barrier (green lithium heparin) tube; SST gel barrier tube; Red top tube Temperature:Refrigerate Place Of Service:PMG-PMP	84403