

LAB	CODE	TEST NAME	CPT4 CODE
PRL	MISC (3117-08)	ABO Genotyping Place of Service:Bloodworks Northwest	
PRL	LAB895	ABO RH Synonyms:ABORH; Blood Type Container:6.0 mL Pink Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OSS, OSV, OWF	
PRL	LAB2244	ABO RH, NEONATAL Container:6.0 mL Pink Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF, PRL BB	
PRL	LAB43	ACETAMINOPHEN LEVEL Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	80307
PRL	LAB836	ACETYLCHOLINE RECEPTOR AB, BINDING Synonyms:Acetylcholine Receptor Binding Antibody; AChR-Binding Antibodies; ACRAb; Myasthenia Gravis Antibody Place of Service:LABCORP BURLINGTON	
PRL	LAB838	ACETYLCHOLINE RECEPTOR AB, BLOCKING Synonyms:AChR-Blocking Antibodies; ACRAb; Anti-AChR Antibody; Myasthenia Gravis Antibody Place of Service:LABCORP BURLINGTON	
PRL	LAB837	ACETYLCHOLINE RECEPTOR AB, MODULATING Place of Service:Labcorp	
PRL	LAB24413	ADAMTS13 Activity Synonyms:0030056; ADAMTS-13; ADAMTS13; MTS13; MTS13 activity Place of Service:ARUP	
PRL	LAB23404	ADAMTS13 Reflex Panel Alert:CRITICAL FROZEN. Synonyms:3000239; ADAMTS-13; ADAMTS13; MTS13; MTS13 activity Container:Light Blue (Sodium Citrate). Refer to Specimen Handling at aruplab.com for hemostasis/thrombosis specimen handling guidelines. Place of Service:ARUP	
PRL	LAB12869	ADIPONECTIN Synonyms:Total Adiponectin Place of Service:LABCORP BURLINGTON	
PRL	LAB511	ADRENOCORTICOTROPIC HORMONE Place of Service:Labcorp	
PRL	LAB901	Aerobe Identification Synonyms:Culture Identification Only; Identification Only, Bacteria; Organism Identification Container:Actively growing organism, in pure culture on agar slant; Also acceptable: Swab in bacterial transport media Place of Service:LabCorp Halsey Microbiology	87077
PRL	LAB16091	AFB MTB COMPLEX, NAAT	
PRL	LAB24370	AGGRESSIVE B-CELL LYMPHOMA, FISH	
PRL	LAB45	ALBUMIN Container:5.0 mL Gold Top Tube (SST); Also acceptable: 6.0 mL Red Top Tube (Plain); 3.0 mL Light Green Top Tube (Lithium Heparin); 4.0 mL Lavender Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	82040