

LAB	CODE	TEST NAME	CPT4 CODE
PRL	LAB2525	Barbiturates Urine Screen Reflex Place of Service:Labcorp	
PRL	LAB15066	BARTONELLA HENSELAE AB, IGG AND IGM Place of Service:LABCORP BURLINGTON	
PRL	LAB15	BASIC METABOLIC PANEL Synonyms:BMP Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	80048
PRL	LAB25811	B-CELL CD20 EXPRESSION Synonyms:CD20 on B Lymphocytes; CD20; % CD20 B Lymphocytes Place of Service:LabCorp Halsey	
PRL	LAB2416	BCR/ABL, P210 QUANT MONITOR Synonyms:BCR-ABL1 Transcript Detection by <RT-PCR, Quantitative; Chronic Myelogenous Leukemia, BCR-ABL1 Transcript Detection by RT-PCR, Quantitative; Quantitative BCR-ABL1 Translocation Detection by RT-PCR for CML Place of Service:LABCORP RTP	
PRL	LAB15032	BCR-ABL BY FISH Synonyms:Philadelphia Translocation; Philadelphia Chromosome; BCR-ABL FISH; BCR/ABL1 t(9; 22); Disease(s): CML, ALL, MPN; Probes: ABL1 (9q34); ASS1 (9q34); BCR (22q11.2); BCRABLFISH (NEO); NeoGenomics; Epic Name: BCR-ABL by Fish Container:6.0 mL Dark Green Top Tube (Sodium Heparin) ; Also acceptable: 4.0 mL Lavender Top Tube (EDTA); Non-diluted bone marrow aspirate collected in a heparinized syringe Place of Service:NeoGenomics Laboratories	88374 automated. Codes may differ if manual analysis is performed
PRL	LAB15039	BCR-ABL1 FOR CML AND ALL Synonyms:Acute Lymphocytic Leukemia, BCR-ABL1 Transcript Detection by RT-PCR, Quantitative; Chronic Myelogenous Leukemia, BCR-ABL1 Transcript Detection by RT-PCR, Quantitative; Quantitative BCR-ABL1 Translocation Detection by RT-PCR for CML and ALL Place of Service:LABCORP RTP	
PRL	LAB24786	BCR-ABL1 KINASE DOMAIN MUTATION ANALYSIS	
PRL	LAB24664	Benzodiazepine Urine Screen Reflex Place of Service:LabCorp	
PRL	LAB24167	BENZODIAZEPINES CONFIRMATION, QUANT, BLOOD Place of Service:LABCORP MEDTOX LABORATORIES INC	
PRL	LAB1831	BETA HYDROXYBUTYRATE, QUANT Container:3.0 mL Light Green Top Tube (Lithium Heparin); Also acceptable: 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain); 4.0 mL Gray Top Tube (Sodium Fluoride/Potassium Oxalate); 4.0 mL Lavender Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	82010
PRL	LAB1179	Beta-2 Glycoprotein 1 Antibodies, IgG and IgM	
PRL	LAB2365	Beta-2 Microglobulin, Serum Synonyms:Beta-2 Microglobulin, Serum Container:5.0 mL Gold Top Tube (SST); Also acceptable: 6.0 mL Red Top Tube (Plain) Place of Service:LabCorp Halsey	82232
PRL	LAB13145	BILE ACIDS, FRACTIONATED Place of Service:LABCORP ESOTERIX ENDOCRINOLOGY	
PRL	LAB693	BILE ACIDS, TOTAL Synonyms:Glyco and Taurochenodeoxycholic Acid Place of Service:LABCORP SEATTLE	