

| LAB | CODE     | TEST NAME   | CPT4 CODE |
|-----|----------|---|-----------|
| PRL | LAB3104  | Lacosamide, Serum<br>Synonyms:Vimpat(R)   |           |
| PRL | LAB96    | LACTATE DEHYDROGENASE<br>Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain)<br>Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV   | 83615     |
| PRL | LAB188   | LACTATE DEHYDROGENASE, BODY FLUID<br>Synonyms:Lactate Dehydrogenase, Variable; LD Variable; LD, Body Fluid  |           |
| PRL | LAB97    | LACTATE DEHYDROGENASE, ISOENZYME, SERUM<br>Synonyms:Lactic Acid Dehydrogenase Isoenzymes; LD Isoenzymes; LDH Isoenzymes<br>Place of Service:LABCORP BURLINGTON  |           |
| PRL | LAB26182 | LACTATE RFLX 2HR IF INDICATED<br>Container:4.0 mL Gray Top Tube (Sodium Fluoride/Potassium Oxalate)<br>Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF  |           |
| PRL | LAB95    | LACTIC ACID<br>Container:4.0 mL Gray Top Tube (Sodium Fluoride/Potassium Oxalate)<br>Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF  | 83605     |
| PRL | LAB1040  | LACTOFERRIN, FECAL, QUANTITATIVE<br>Synonyms:Fecal Lactoferrin; Fecal Leukocytes; IBD-Scan(R); Stool Lactoferrin<br>Place of Service:LABCORP SPOKANE  |           |
| PRL | LAB475   | Lamotrigine<br>Synonyms:Lamictal(R)   |           |
| PRL | LAB830   | LEAD, CAPILLARY BLOOD<br>Synonyms:Pb, Blood<br>Place of Service:LABCORP SPOKANE   |           |
| PRL | LAB12997 | LEAD, INDUSTRIAL EXPOSURE PANEL<br>Synonyms:Lead and Protoporphyrin (ZPP)<br>Place of Service:LABCORP SPOKANE   |           |
| PRL | LAB98    | LEAD, WHOLE BLOOD<br>Synonyms:Pb, Blood<br>Place of Service:LABCORP SPOKANE   |           |
| PRL | LAB886   | Legionella Antigen, Qualitative<br>Synonyms:Legionella Urinary Antigen<br>Container:Sterile, Leak-Proof Container ; Also acceptable: BD Vacutainer® Urine Boric Acid Tube<br>Place of Service:LabCorp Halsey  | 87899     |
| PRL | LAB13228 | Legionella, NAAT<br>Synonyms:Legionella PCR; Legionella species PCR; LEG PCR<br>Container:BD Universal Viral Transport System (UTM); Also acceptable: Viral Transport Media (VTM); Sterile Leak-Proof Container<br>Place of Service:LabCorp Halsey Microbiology | 87541     |
| PRL | LAB13229 | LEPTIN<br>Place of Service:LABCORP BURLINGTON   |           |
| PRL | LAB99    | LIPASE<br>Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain)<br>Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF   | 83690     |