

Online Test Menu **Oregon Regional Laboratories**



		Oregon Regional Eaboratories	•
LAB	CODE	TEST NAME	CPT4 CODE
PRL	LAB103	MAGNESIUM Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	83735
PRL	LAB13238	MAGNESIUM, PERINATAL Container:5.0mL Gold Top Tube (SST); 6.0 mL Plain Red Top, 3.0 mL Light Green Top Tube (PST) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV	
PRL	LAB2404	MAGNESIUM, RBC Synonyms:Mg, RBC Place of Service:LABCORP BURLINGTON	
PRL	LAB406	MAGNESIUM, URINE, 24HR Place of Service:Labcorp	
PRL	LAB23147	MALARIA SCREEN WITH REFLEX TO QUANT AND SPECIATION, BLOOD Alert:Send STAT to lab due to short stability and due to critical nature of the test Sample must be RECEIVED in lab within 24 hours of collection Container:4.0 mL Lavender Top Tube (EDTA) Also acceptable: 6.0 mL Pink Top Tube (EDTA) Capillary blood obtained by finger stick with smears prepared at bedside Place of Service:OHR, OMW, ONB, OPH, OSS, OSV, OWF, Asante Rogue RMCL (OMD)	87207 87015 87899 87207-26
PRL	LAB883	MALARIA, SMEAR Container:4.0 mL Lavender Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	
PRL	LAB1050	MANGANESE Synonyms:Mn, Plasma Place of Service:NMS via LabCorp	
PRL	LAB1052	MANGANESE, WHOLE BLOOD Synonyms:Mn, Blood Place of Service:LABCORP SPOKANE	
PRL	LAB657	Measles (Rubeola) Antibody, IgG Synonyms:Measles	
PRL	LAB408	MERCURY, URINE, 24HR Synonyms:Hg, Urine Place of Service:LABCORP SPOKANE	
PRL	LAB831	MERCURY, WHOLE BLOOD Synonyms:Hg, Blood Place of Service:LABCORP SPOKANE	
PRL	LAB2135	METANEPHRINES, FRACTIONATED, FREE, PLASMA Synonyms:Metanephrine; Normetanephrine; Plasma Metanephrines; Quantitative Metanephrines Place of Service:LABCORP SPOKANE	
PRL	LAB12316	METANEPHRINES, FRACTIONATED, URINE, 24HR Place of Service:LABCORP SPOKANE	
PRL	LAB409	METANEPHRINES, TOTAL, URINE Place of Service:LABCORP ESOTERIX ENDOCRINOLOGY	
PRL	LAB15018	Methanol, Blood Alert:For OMD, OSS, OHR ONLY! All other sites use: MISC (007062) to Legacy via LabCorp.	