

LAB	CODE	TEST NAME	CPT4 CODE
PRL	LAB748	UREA NITROGEN, URINE Container:Random Urine in a sterile, leak-proof container (urine cup) without preservative Place of Service:LabCorp Halsey	84540
PRL	LAB123748	Urea Nitrogen, Urine Timed Synonyms:BUN; Nitrogen Balance Study; Urine Urea Nitrogen Container:24-hour or timed urine in a sterile leak-proof urine container without preservative Place of Service:LabCorp Halsey	84540
PRL	LAB23927	UREAPLASMA/MYCOPLASMA HOMINIS Synonyms:Culture, Mycoplasma hominis (genital specimens); Culture, Ureaplasma urealyticum, Genital; Mycoplasma hominis Genital Culture; Mycoplasma T-Strain Culture, Genital; Ureaplasma Culture Place of Service:LABCORP BURLINGTON	
PRL	LAB141	URIC ACID Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain); 4.0 mL Lavender Top Tube (EDTA) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	84550
PRL	LAB11450	Uric Acid, Urine Random Container:Random Urine in a sterile, leak-proof container (urine cup) without preservative Place of Service:LabCorp Halsey	84560
PRL	LAB841	Uric Acid, Urine Timed Container:24-hour or timed urine in a sterile leak-proof urine container without preservative Place of Service:LabCorp Halsey	84560
PRL	LAB347	URINALYSIS Container:Clean Catch Urine in a sterile, leak-proof container (urine cup) without preservative; Also acceptable: Urinalysis Preservative Tube (Marble Top Tube) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF, OPH PCI NEWBERG, OWF PCI WILLAMETTE FALLS	81003
PRL	LAB348	URINALYSIS WITH MICROSCOPIC Container:Clean Catch Urine in a sterile, leak-proof container (urine cup) without preservative; Also acceptable: Urinalysis Preservative Tube (Marble Top Tube) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	81001
PRL	LAB2479	URINALYSIS WITH MICROSCOPIC IF INDICATED Container:Clean Catch Urine in a sterile, leak-proof container (urine cup) without preservative; Also acceptable: Urinalysis Preservative Tube (Marble Top Tube) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	
PRL	LAB2209	URINALYSIS, MICROSCOPIC ONLY Container:Clean Catch Urine in a sterile, leak-proof container (urine cup) without preservative; Also acceptable: Urinalysis Preservative Tube (Marble Top Tube) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	81015
PRL	LAB12834	URINALYSIS, REFLEX MICROSCOPIC AND/OR CULTURE Container:Clean Catch Urine in a sterile, leak-proof container (urine cup) without preservative; Also acceptable: Urinalysis Preservative Tube (Marble Top Tube) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	81003 81001 if microscopic added. Culture and susceptibility tests are billed separately from Urinalysis