

Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SED	LAB262	NEISSERIA GONORRHOEAE, NAAT Container:BD Probe Tec Swab Collection Kit for cervical and urethral specimens or a Screw-cap sterile cup for urine (15-20 mL) specimens	87591
SED	LAB892	NEONATAL BLOOD WORKUP Container:Pink top K2EDTA (6 mL)	86900 86901 86880
SED	LAB23918	NEUROMYELITIS OPTICA AB, SERUM Container:SST gel-barrier tube Temperature:Refrigerate	86053
SED	LAB13263	NEURON SPECIFIC ENOLASE AG	
SED	LAB13262	NEURON SPECIFIC ENOLASE, CSF	
SED	LAB2150	NEUTROPHIL OXIDATIVE BURST Container:Heparinized Green Top Tube (1 mL) Temperature:Ambient	88184 88185 88187
SED	LAB480	NEWBORN BLOOD SPOT SCREENING Container:Newborn Screening Collection Cards	multiple
SED	LAB1795	NEWBORN METABOLIC SCREEN REPEAT Container:Newborn Screening Collection Cards	no charges
SED	LAB2435	NICOTINE AND METABOLITES	
SED	LAB414	NICOTINE, URINE, QUAL	
SED	LAB24278	NMR LIPOPROFILE TEST (LP-IR VALUES)	
SED	LAB13270	NOROVIRUS 1 AND 2, NAAT	
SED	LAB816	N-TELOPEPTIDE, CROSS-LINKED, URINE	
SED	LAB1780	NT-PRO BNP	
LABCORP	183624	N. gonorrhoeae NAA, Confirm Container:Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) or SurePath(TM) liquid cytology vials Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87591