

LAB	CODE	TEST NAME	CPT4 CODE
PSMMC	LAB24382	RAINBOW BLOOD PANEL Container:Blue (Sodium Citrate) tube, Gold SST, 2 Green (Lithium Heparin), Lavender (EDTA), Grey (Na Fluoride) tube Place Of Service:PSMMC	
PSMMC	LAB23143	Rapid Plasma Reagin (RPR), Quantitation Place Of Service:LABCORP	
PSMMC	LAB494	RAPID PLASMA REAGIN, QUAL Container:Gold SST or Red-top tube Temperature:Refrigerate Place Of Service:PSMMC	86592
PSMMC	LAB19	RENAL FUNCTION PANEL Container:Green (Lithium Heparin) tube, Gold SST or Red-top tube Place Of Service:PSMMC	80069
PSMMC	LAB532	RENIN ACTIVITY Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB1307	Respiratory Pathogen Panel, NAAT Container:Place swab into a vial containing viral transport medium ("VTM": e.g. M4, M4-RT, M5, M6; Universal Transport Media; and Universal Viral Transport Media). Temperature:Specimens used for RP Flex testing must be tested or stored at 2-8°C within 4 hours of collection. Place Of Service:PSMMC	87633 87798 x3
PSMMC	LAB296	RETIC COUNT Container:Lavender EDTA Tube Temperature:Refrigerate Place Of Service:PSMMC	85045
PSMMC	LAB773	RETICULIN AB Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB23639	RETINOL BINDING PROTEIN Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB23684	RHEUMATOID FACTOR, IGA, IGG AND IGM Place Of Service:ESOTERIX ENDOCRINOLOGY CF	
PSMMC	LAB206	RHEUMATOID FACTOR, QUANT Place Of Service:LABCORP SEATTLE	
PSMMC	LAB1399	RIBOSOMAL P AB, IGG Place Of Service:LABCORP SPOKANE	
PSMMC	LAB1264	RICKETTSIA RICKETTSII AB Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB13292	RNA POLYMERASE III AB IGG Place Of Service:ESOTERIX ENDOCRINOLOGY CF	
PSMMC	LAB772	RNP AB, IGG Place Of Service:LABCORP SPOKANE	