

LAB	CODE	TEST NAME	CPT4 CODE
PSMMC	LAB5687	VAGINAL PATHOGENS DNA DIRECT PROBE Container:Affirm VPIII Ambient Temperature Transport System (ATTS) Temperature:Maintain specimen at room temperature. Place Of Service:PSMMC	87480 87510 87660
PSMMC	LAB10042	VAGINITIS WITH CANDIDA Place Of Service:InCyte	
PSMMC	LAB24	VALPROIC ACID LEVEL Container:Green (Lithium Heparin) tube, Gold SST or Red-top tube Temperature:Refrigerate Place Of Service:PSMMC	80164
PSMMC	LAB40	VANCOMYCIN LEVEL Container:Green (Lithium Heparin) tube, Gold SST or Red-top tube Temperature:Refrigerate Place Of Service:PSMMC	80202
PSMMC	LAB41	VANCOMYCIN, PEAK Container:Green (Lithium Heparin) tube, Gold SST or Red-top tube Temperature:Refrigerate Place Of Service:PSMMC	80202
PSMMC	LAB39	VANCOMYCIN, TROUGH Container:Green (Lithium Heparin) tube, Gold SST or Red-top tube Temperature:Refrigerate Place Of Service:PSMMC	80202
PSMMC	LAB452	VANILLYLMANDELIC ACID, URINE Place Of Service:LABCORP	
PSMMC	LAB750	VANILLYLMANDELIC ACID, URINE, RANDOM Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB162	VARICELLA ZOSTER AB, IGG Place Of Service:LABCORP SEATTLE	
PSMMC	LAB22164	VARICELLA ZOSTER AB, IGG AND IGM Place Of Service:LABCORP SPOKANE	
PSMMC	LAB163	VARICELLA ZOSTER AB, IGM Place Of Service:LABCORP SEATTLE	
PSMMC	LAB2177	VARICELLA ZOSTER, NAAT Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB24765	VASCULAR ENDOTHELIAL GROWTH FACTOR Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB1107	VASOACTIVE INTESTINAL PEPTIDE Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB859	VDRL, CSF Place Of Service:LABCORP BURLINGTON	