

LAB	CODE	TEST NAME	CPT4 CODE
PSMMC	LAB1271	YERSINIA ANTIBODY, IGA, IGM, IGG Place Of Service:ARUP LABORATORIES INC	
LABCORP	512116	Y Chromosome Microdeletion Analysis Container:Lavender-top (EDTA) tube or pink-top (EDTA) tube or tan-top (EDTA) tube or yellow-top (ACD) tube or PurFlock buccal swab kit or Oragene Dx 500 saliva collection kit. Temperature:Maintain specimen at room temperature or refrigerate at 4C. Do NOT freeze.	81403
LABCORP	183255	Yeast Suscep, Isavuconazole Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	182329	Yeast Susceptibility Testing, Amphotericin Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183483	Yeast Susceptibility Testing, Anidulafungin Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183491	Yeast Susceptibility Testing, Caspofungin Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183257	Yeast Susceptibility Testing, Comprehensive Panel Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	182220	Yeast Susceptibility Testing, Fluconazole Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	182436	Yeast Susceptibility Testing, Itraconazole Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183509	Yeast Susceptibility Testing, Micafungin Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183517	Yeast Susceptibility Testing, Posaconazole Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183260	Yeast Susceptibility Testing, Rezafungin Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	183525	Yeast Susceptibility Testing, Voriconazole Container:A swab transport inoculated with the isolate from a pure culture or an isolate growing on an agar slant in screw-cap container packed as an etiologic agent. Maintain a subculture at the submitting laboratory. Temperature:Maintain media at room temperature.	87186
LABCORP	829957	Yersinia IgM by Immunoblot Container:Red-top tube or gel-barrier tube	86793

LAB	CODE	TEST NAME	CPT4 CODE
-----	------	-----------	-----------

Temperature:Refrigerate