

LAB	CODE	TEST NAME	CPT4 CODE
PMC	LAB720	17-OH Progesterone LCMS	
PMC	LAB410	24 Hour Microalbumin Creatinine Ratio	
PMC	LAB9089	24 Hour Urine Calcium Temperature:keep refrigerated	
PMC	LAB9198	24 Hour Urine Magnesium Temperature:keep refrigerated	
PMC	LAB841	24 Hour Urine Uric Acid Temperature:keep refrigerated	
PMC	LAB10080	25-Hydroxyvitamin D LCMS D2+D3	
PMC	LAB14095	4-Plex Sars CoV-2, Flu A, Flu B, RSV Container:VTM/UTM	
PMC	LAB352	5-HIAA,Quant.,24 Hr Urine	
LABCORP	071910	071910 7 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	074666	074666 9 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	500600	1,25-Dihydroxy,Vitamin D by MS Container:Gel-barrier tube, lavender-top (EDTA) tube, or green-top (heparin) tube Temperature:FREEZE Stability: Room temperature up to 2 days Refrigerated up to 2 days Frozen up to 544 days Freeze/Thaw cycle = 6	82652
LABCORP	791000	1,3-Dimethylamylamine, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for eight days at room temperature, refrigerated, or frozen.	80307
LABCORP	761490	10+Oxycodone-Scr, Oral Fluid Container:Quantisal(TM) Oral Fluid collection device, contact lab for collection kit. Temperature:Maintain specimen at room temperature for transport to laboratory; refrigerate if transport is delayed.	80307
LABCORP	768561	11+Oxyco+Crt-Bund+1-Confirm Container:Urine container Temperature:Maintain specimen at room temperature Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80346 80307
LABCORP	500171	11-Deoxycortisol: HPLC/MS-MS Container:Gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube Temperature:Freeze	82634