

LAB	CODE	TEST NAME	CPT4 CODE
PMC	LAB571	G6PD,Qn,Bld and Red Cell Count	
PMC	LAB650	GAD-65 Autoantibody	
PMC	LAB85	Gamma Glutmyl Transferase (GGT) Container:Lithium Heparin, SST	
PMC	LAB696	Gastric Occult Blood	
PMC	LAB69601	Gastric WBC's	
PMC	LAB80	Gastrin, Serum	
PMC	LAB997101	GC/CT by PCR Container:Voided urine, swabs from endocervical, vaginal, rectal, and pharyngeal sources	
PMC	LAB9561	Genital Culture Container:Swab transported in Culturette	
PMC	LAB14615	Genital Mycoplasma Profile, NAA	
PMC	LAB9460	GenoSure(R) MG	
PMC	LAB28	Gentamicin, Peak Level Container:Plain red top tube or green top tube, or purple top tube, no gel. Temperature:2-8°C. Also acceptable: plasma (heparin) or EDTA. Submit specimen in a Standard Transport Tube.	
PMC	LAB27	Gentamicin, Random Container:Plain red top tube, purple top tube or green top tube no gel. Temperature:2-8°C.. Also acceptable: EDTA or heparin plasma . Submit specimen in a Standard Transport Tube.	80170
PMC	LAB26	Gentamicin, Trough Level Container:Plain red top or purple top tube or green top tube, no gel. Temperature:2-8°C. Also acceptable: plasma (heparin) or EDTA (purple). Submit specimen in a Standard Transport Tube.	80170
PMC	LAB258	Giardia Antigen by EIA Container:Preserve in 10% formalin within one hour of collection. Temperature:Transport stool preserved in 10% formalin at 20-25°C. (Min: 1 g) Submit specimen according to Biological Substance, Category B, shipping guidelines	87329
PMC	LAB11747	Gliadin IgG/IgA Ab Prof, EIA	