

LAB	CODE	TEST NAME	CPT4 CODE
PMC	LAB886	L. pneumophila Serogp 1 Ur Ag	
PMC	LAB95	Lactate Container:One 4 ml tube Lithium heparin,no gel. (Min: 3 mL) Collect on ice at the bedside or colleciton site Draw the specimen without the use of a tourniquet or within three minutes of applying the tourniquet, but before releasing the tourniquet. Temperature:2-8°C.	83605
PMC	LAB97	Lactate Dehydrogenase, Isoenzymes Container:Serum separator tube. Temperature:Room temperature.	83625 83615
PMC	LAB96	Lactate Dehydrogenase, Serum or Plasma Container:One 4 mL plasma separator tube or serum separator tube. Temperature:2-8°C	83615
PMC	LAB731	Lactoferrin, Fecal, Quant.	
PMC	LAB475	Lamotrigine (Lamictal), Serum	
PMC	LAB188	LD, Body Fluid	
PMC	LAB9807	Lead Standard Profile, Blood	
PMC	LAB98	Lead, Blood (Adult)	
PMC	LAB9821	Leflunomide	
PMC	LAB12663	Legionella RT PCR Panel (BAL)	
PMC	LAB902	Legionella Species Culture	
PMC	LAB477	Levetiracetam (Keppra), S	
PMC	LAB99	Lipase, Serum or Plasma Container:One 4 mL plasma separator tube or serum separator tube. Temperature:2-8°C. Submit specimen in a Standard Transport Tube.	83690
PMC	LAB18	Lipid Panel Container:One 4 mL plasma separator tube or serum separator tube. Patient must fast 12-15 hours prior to collection. Temperature:2-8°C. Submit specimen in a Standard Transport Tube.	80061