

LAB	CODE	TEST NAME	CPT4 CODE
SAH	ABRH2	ABO And RH Container:Pink top (PNK) Temperature:Specimens should be tested as soon as possible. Samples that cannot be tested immediately should be stored at 2 degrees to 8 degrees Celsius. Do NOT freeze red cells. Place Of Service:MLAB	86900
SAH	ACETM	ACETAMINOPHEN Container:Mint green top (lithium heparin)(LGR) Temperature:Room Temperature Place Of Service:MLAB	80329
SAH	ACETY	ACHR BINDING ABS Place Of Service:LabCorp	
SAH	ACP	ACID PHOSPTASE Container:Gold or Tiger top(SST) Temperature:FREEZE up to 1 month Place Of Service:MLAB	84066
SAH	ACTH	ACTH Place Of Service:LabCorp	
SAH	ASMA	ACTIN (SMMUSCLE) AB Place Of Service:LabCorp	
SAH	ADA13	ADAMTS13 ACTIVITY Place Of Service:LabCorp	
SAH	AFCXS	AFB CULTURE&SMEAR Place Of Service:LabCorp	
SAH	AFBLD	AFB CULTURE&SMEAR(BLOOD) Place Of Service:LabCorp	
SAH	AFBSF	AFB CULTURE&SMEAR(CSF) Place Of Service:LabCorp	
SAH	AFPT	AFP TETRA Place Of Service:LabCorp	
SAH	AFP	AFP TUMR MKR Place Of Service:LabCorp	
SAH	ALB	ALBUMIN Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	82040
SAH	FLALB	ALBUMIN FLUID Container:Sterile container Temperature:Refrigerate Place Of Service:MLAB	82040
SAH	ALC	ALCOHOL Container:Mint green top (lithium heparin)(LGR) Temperature:Room Temperature Place Of Service:MLAB	80320