

LAB	CODE	TEST NAME	CPT4 CODE
SAH	HPSAG	H. PYLORI STOOL AG Place Of Service:LabCorp	
SAH	HALDOL	HALDOL/HALDOLPERIDOL Place Of Service:LabCorp	
SAH	HAPT	HAPTOGLOBIN Place Of Service:LabCorp	
SAH	HBVQN	HBV REAL-TIME PCR, QUANT Place Of Service:LabCorp	
SAH	HCGQ	HCG QUANT Container:Mint green top (lithium heparin)(LGR) Temperature:Room temperature Place Of Service:MLAB	84702
SAH	RIBA	HCV AB VERIFICATION Place Of Service:LabCorp	
SAH	HCPCR	HCV RT-PCR, QUANT Place Of Service:LabCorp	
SAH	HDL	HDL CHOLESTEROL Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature.[1] Place Of Service:MLAB	83718
SAH	HEEL	Heelstick ABORH Container:Lavender top (EDTA) Micro Tube Place Of Service:MLAB	86900
SAH	HTLYM	HELPER T-LYMPH-CD4 Place Of Service:LabCorp	
SAH	HCT	HEMATOCRIT Container:Lavender top(LAV) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	85014
SAH	HGB	HEMOGLOBIN Container:Lavender top(LAV) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	85018
SAH	HGBEP	HEMOGLOBINOPATHY PROFILE Place Of Service:LabCorp	
SAH	HEMOG	HEMOGRAM Container:Lavender top(LAV) Temperature:Maintain specimen at room temperature up to 24 hours. Place Of Service:MLAB	85025
SAH	HBAAT	HEP A AB TOTAL Place Of Service:LabCorp	