

LAB	CODE	TEST NAME	CPT4 CODE
SAH	NICOT	NICOTINE AND METABOLITE, UR QN Place Of Service:LabCorp	
SAH	NCOMP	NOCARDIA/ACTNO COMPLEX Place Of Service:LabCorp	
LABCORP	183624	N. gonorrhoeae NAA, Confirm Container:Gen-Probe(R) APTIMA(R) swab or APTIMA(R) urine specimen transport; ThinPrep(R) or SurePath(TM) liquid cytology vials Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87591
LABCORP	247315	N.Meningitidis A/C/W-135/Y IgG Container:Gel-barrier tube or red-top tube Temperature:Refrigerated	86741 x4
LABCORP	222000	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84300 84133 82436
LABCORP	205914	Na U+Cl U+K U Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84300 84133 82436
LABCORP	810960	Nabumetone as Metabolite Container:Red-top tube, lavender-top (EDTA) tube fluoride) tube Temperature:Refrigerated Stability: Room Temperature: 30 days Refrigerated: 30 days Frozen: 17 months	80299
LABCORP	831560	N-Acetyl-L-Aspartic Acid,Urine Container:Plastic urine container Temperature:FREEZE	82542
LABCORP	805241	NailStat-12 Container:Contact lab for Hair collection kit. Temperature:Maintain specimen at room temperature.	80307
LABCORP	805131	NailStat-14 Container:Contact lab for Hair collection kit. Temperature:Maintain specimen at room temperature.	80307
LABCORP	805239	NailStat-9 Container:Contact lab for Hair collection kit. Temperature:Maintain specimen at room temperature.	80307
LABCORP	764175	Nalbuphine Container:Urine container Temperature:Maintain specimen at room temperature Stability: If arrival will extend beyond 7 days then refrigerate.	80307
LABCORP	701195	NALOXONE,UR,MS CONFIRM Container:Plastic urine container without preservative Temperature:Room Temperature	80362
LABCORP	702655	Naltrexone and Mtb, MS, UR Container:Urine Bottle	80362
LABCORP	702518	Naltrexone, MS Confirm, Urine, Medication Assisted Treatment(MAT) Container:Plastic urine container without preservative. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80362