

LAB	MNEMONIC	TEST NAME	CPT CODE
SPH	PANCPOLY	Pancreatic Polypeptide	
SPH	PAPHPVAG	PAP W/HPV	
SPH	PAPIGHPVCN	PAP,CTNG W/HPV	
SPH	PAPAGEGD	PAP,CTNG W/HPV (AGE/RISK BASED)	
SPH	PAPIGHPCNT	PAP,CTNGTV&HPV(AGE/RISK Based)	
SPH	PTT	Partial Thromboplastin Time Container:Blue (Sodium Citrate) Temperature:Room Temperature: up to 4 hours from time of collection Frozen: remove plasma from cells. Freeze plasma for up to two weeks	85730
SPH	PARVOB19	Parvovirus B19, Human, IgG/IgM	
SPH	PATH SMEAR REV	Path Smear Review (SPH) Container:Viral Transport Media	
SPH	PATHEXAM3	Pathology Examination Container:Sterile Cup Temperature:Room Temperature	
SPH	PTFPNL	Peritoneal Fluid Panel Container:Syringe Temperature:Room Temperature	89051
SPH	PERT PCR	Pertussis Bordetella PCR	
SPH	UPH	pH Urine Container:Urinalysis Tube Temperature:Refrigerated	81003
SPH	PHENO	Phenobarbital, Serum	
SPH	PHENYTOIN	Phenytoin Level (Reference) Container:Green (Lithium Heparin) Temperature:See SSM test catalog	80185
SPH	PHOSIGA	Phospholipid Ab, IgA	