

Online Test Menu Laboratory Catalog



| LAB | CODE | TEST NAME | CPT4 CODE |
|-----|----------|---|-----------|
| SPH | LAB13117 | ABO ONLY Container:Pink (preferred); Lavender, Purple, Red Temperature:Refrigerated 2-8° C | |
| SPH | LAB895 | ABO RH Container:Pink (preferred); Lavender, Purple, Red; Must use blood band if transfusing patient Temperature:Refrigerated 2-8° C | |
| SPH | LAB2244 | ABO RH, NEONATAL Container:Pink (preferred); Lavender, Purple, Red; Must use blood band if transfusing patient Temperature:Refrigerated 2-8° C | |
| SPH | LAB43 | ACETAMINOPHEN LEVEL Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab. | 80307 |
| SPH | LAB456 | ACETONE Place Of Service:LabCorp | |
| SPH | LAB836 | ACETYLCHOLINE RECEPTOR AB, BINDING Place Of Service:LabCorp | |
| SPH | LAB838 | ACETYLCHOLINE RECEPTOR AB, BLOCKING Place Of Service:LabCorp | |
| SPH | LAB428 | ACETYLCHOLINESTERASE Place Of Service:LabCorp | |
| SPH | LAB429 | ACYLCARNITINES, PLASMA, QUANT Place Of Service:LabCorp | |
| SPH | LAB24699 | ADALIMUMAB AND ANTIBODIES Place Of Service:LabCorp | |
| SPH | LAB23404 | ADAMTS 13 ACTIVITY , REFLEX Place Of Service:LabCorp | |
| SPH | LAB24413 | ADAMTS13 EVALUATION Place Of Service:LabCorp | |
| SPH | LAB24290 | ADENOSINE DEAMINASE, PLEURAL Place Of Service:LabCorp | |
| SPH | LAB2530 | ADENOVIRUS DNA, QUANTITATIVE, RT-PCR Place Of Service:VIRACOR EUROFINS | |
| SPH | LAB511 | ADRENOCORTICOTROPIC HORMONE Place Of Service:LABCORP SEATTLE | |