

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB24734	N METHYL D ASPARTATE RECPT IGG Place Of Service:LabCorp	
SPH	LAB13255	NATURAL KILLER CELL Place Of Service:LabCorp	
SPH	LAB262	NEISSERIA GONORRHOEAE, NAAT Container:Endocervical swab specimens, patient-collected vaginal swab specimens, oropharyngeal (throat) swabs, and rectal swabs must be collected using the Xpert Swab Specimen Collection Kit, first-catch urine may be collected in a sterile container or Xpert CT/NG Urine Specimen Collection Kit Temperature:room temperature or refrigerated Place Of Service:Testing platform is Cepheid GeneXpert. Acceptable specimens are vaginal and endocervical swabs, oropharyngeal (throat) swabs, rectal swabs, and first-catch urine. Please call lab prior to collection for instructions and collection materials for reference lab testing on urethral swabs or patients under the age of 14, which can be performed via the Aptima platform.	
SPH	LAB2193	NEUROMYELITIS OPTICA AB Place Of Service:LABCORP BURLINGTON	
SPH	LAB13263	NEURON SPECIFIC ENOLASE AG Place Of Service:LabCorp	
SPH	LAB23621	NEURONAL CELL AB Place Of Service:LabCorp	
SPH	LAB1795	NEWBORN METABOLIC SCREEN REPEAT Place Of Service:LabCorp	
SPH	LAB2435	NICOTINE AND METABOLITES Place Of Service:LabCorp	
SPH	LAB24278	NMR LIPOPROFILE TEST (LP-IR VALUES) Place Of Service:LabCorp	
SPH	LAB13270	NOROVIRUS 1 AND 2, NAAT Place Of Service:LabCorp	
SPH	LAB807	NORTRIPTYLINE LEVEL Place Of Service:LabCorp	
SPH	LAB816	N-TELOPEPTIDE, CROSS-LINKED, URINE Place Of Service:LabCorp	
SPH	LAB1780	NT-PRO BNP Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	83880
SPH	LAB23624	NUCLEOPHOSMIN (NPM1) MUTATION ANALYSIS PCR AND FRAGMENT ANALYSIS Place Of Service:GenPath	
LABCORP	831560	N-Acetyl-L-Aspartic Acid,Urine Container:Plastic urine container Temperature:FREEZE	82542