

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB494	RAPID PLASMA REAGIN, QUAL Place Of Service:LabCorp	
SPH	LAB19	RENAL FUNCTION PANEL Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	80069
SPH	LAB532	RENIN ACTIVITY Place Of Service:LABCORP BURLINGTON	
SPH	LAB1861	RESPIRATORY VIRUS PANEL NAAT Container:VTM": e.g. M4, M4-RT, M5, M6; Universal Transport Media; and Universal Viral Transport Media Temperature:Room temperature (18-30 C) for up to 4 hours, at refrigerator temperature (2-8 C) for up to 3 days, or at freezer temperature (<15 C) for up to 30 days.	87998 x3 87633 x3
SPH	LAB296	RETIC COUNT Container:Purple Top	
SPH	LAB773	RETICULIN AB Place Of Service:LabCorp	
SPH	LAB23639	RETINOL BINDING PROTEIN Place Of Service:LabCorp	
SPH	LAB1787	RH ONLY Container:Pink top preferred; Lavender, Purple, Red Temperature:Refrigerated 2-8° C	
SPH	LAB23684	RHEUMATOID FACTOR, IGA, IGG AND IGM Place Of Service:LabCorp	
SPH	LAB206	RHEUMATOID FACTOR, QUANT Place Of Service:LabCorp	
SPH	LAB1399	RIBOSOMAL P AB, IGG Place Of Service:LabCorp	
SPH	LAB1264	RICKETTSIA RICKETTSII AB Place Of Service:LabCorp	
SPH	LAB13292	RNA POLYMERASE III AB IGG Place Of Service:LabCorp	
SPH	LAB772	RNP AB, IGG Place Of Service:LabCorp	
SPH	LAB2321	ROTAVIRUS AG, EIA Place Of Service:LabCorp	