

Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB13043	VAGINAL PATHOGENS DNA Container:Xpert Swab specimen collection kit (pink swab) Temperature:Room Temperature, refrigerated	
SPH	LAB5687	VAGINAL PATHOGENS DNA DIRECT PROBE Place Of Service:LabCorp	
SPH	LAB24	VALPROIC ACID LEVEL Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	80164
SPH	LAB40	VANCOMYCIN LEVEL Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	80202
SPH	LAB41	VANCOMYCIN, PEAK Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	80202
SPH	LAB39	VANCOMYCIN, TROUGH Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	80202
SPH	LAB750	VANILLYLMANDELIC ACID, URINE, RANDOM Place Of Service:LabCorp	
SPH	LAB162	VARICELLA ZOSTER AB, IGG Place Of Service:LabCorp	
SPH	LAB22164	VARICELLA ZOSTER AB, IGG AND IGM Place Of Service:LabCorp	
SPH	LAB163	VARICELLA ZOSTER AB, IGM Place Of Service:LabCorp	
SPH	LAB2177	VARICELLA ZOSTER, NAAT Place Of Service:LabCorp	
SPH	LAB24765	VASCULAR ENDOTHELIAL GROWTH FACTOR Place Of Service:LabCorp	
SPH	LAB1107	VASOACTIVE INTESTINAL PEPTIDE Place Of Service:LabCorp	
SPH	LAB859	VDRL, CSF Place Of Service:LabCorp	
SPH	LAB765	VISCOSITY, SERUM Place Of Service:LABCORP BURLINGTON	