

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
BSHSI-R	ZIKML	ZIKA VIRUS AB, IGM	
BSHSI-R	ZINCL	ZINC	
BSHSI-R	ZINWL	ZINC, WHOLE BLOOD	
BSHSI-R	ZNT8L	ZNT8 AUTO-ANTIBODIES	
BSHSI-R	ZONIL	ZONISAMIDE(ZONEGRAN)	
LABCORP	790660	Zaleplon, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LABCORP	790710	Zaleplon, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80368
LABCORP	010180	Zinc Protoporphyrin (ZPP), Workplace Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube. Temperature:Room temperature. Stable for five days at room temperature or for seven days when refrigerated.	84202
LABCORP	070111	Zinc Random Urine	84630 82570
LABCORP	070029	Zinc, RBC Container:Royal blue-top (EDTA) tube Temperature:Room Temperature	84630
LABCORP	003434	Zinc, Urine Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84630 82570
LABCORP	790384	Ziprasidone, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:Maintain specimen at room temperature; stable for three days. Stable refrigerated 14 days	80299
LABCORP	808135	Zolpidem, Serum or Plasma Container:Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for three days at room temperature, for three days refrigerated, and for two weeks frozen.	80299
LABCORP	081100	Zolpidem, Ur Container:Plastic urine container without preservative Temperature:Maintain specimen at room temperature NOTE: Room temperature is adequate for shipping and/or short-term storage (up to 3 days); specimens may also be refrigerated or frozen	80368
LABCORP	809646	Zolpidem, Whole Blood Container:Lavender Top (EDTA) Tube OR Gray Top (Sodium Fluoride) Tube OR Green Top (Sodium Heparin) Tube. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80368