



Online Test Menu



LABORATORY CATALOG

Where good health begins.

| LAB | CODE | TEST NAME | CPT4 CODE |
|-----|---------|--|-----------|
| LRH | PALIP | PALIPERIDONE | |
| LRH | PANCSTT | PANCREASTATIN | |
| LRH | PANCELS | PANCREATIC ELASTASE-1 | |
| LRH | PANCP0 | PANCREATIC POLYPEPTIDE | |
| LRH | PARASIT | PARASITE EXAMINATION, BLOOD | |
| LRH | WORMID | PARASITE ID, WORM | |
| LRH | PARIETA | PARIETAL CELL TOTAL AUTOABS | |
| LRH | PARVBGM | PARVOVIRUS B19 IGG & IGM ABS | |
| LRH | PARVPCR | PARVOVIRUS B19 QUANT PCR | |
| LRH | PH | PH (BLOOD) Container:Sodium or Lithium heparin syringe or Capillary collection device; Venous blood in green-top heparinized tube. Temperature:Ambient | 82800 |
| LRH | BFPH | PH, BODY FLUID | |
| LRH | FECPH | PH, FECES | |
| LRH | UPH | PH, URINE Container:Chemically clean or sterile disposable specimen cup with the leak proof lid. Pediatric collections are performed using clear polyethylene bags. Temperature:ambient or refrigerate | 83986 |
| LRH | UPPCPF | PHENCYCLIDINE CONFIRM. URINE | G0480 |
| LRH | PHENB | PHENOBARBITAL Container:Light Green (PST) tube; Green-top (heparin) tube, Serum separator-tube, or Red-Top Tube Temperature:Ambient, Refrigerated, or Frozen | 80184 |