

LAB	CODE	TEST NAME	CPT4 CODE
LRH	SACCHA	SACCHAROMYCES CEREVISIAE IGG & IGA	
LRH	SALIC	SALICYLATE Container:Light Green (PST) tube; Green-top (heparin) tube, Serum separator-tube, or Red-Top Tube Temperature:Ambient, Refrigerate, or Frozen	80307
LRH	SARSSF	SARS ANTIGEN SOFIA	
LRH	SARSG2	SARS-CoV-2 IgG II ANTIBODY, LRH	
LRH	SARSGSQ	SARS-CoV-2 Semi-Quantitative Total Antibody, Spike	
LRH	COVID	SARS-CoV-2, NAA	
LRH	SCL70	SCL-70 IgG AUTOABS	
LRH	ESR	SEDIMENTATION RATE Container:Lavender top (EDTA) tube Temperature:Ambient or Refrigerated	85652
LRH	SELEN	SELENIUM	
LRH	SEMENPV	SEMEN - POST VASECTOMY SPERM COUNT Container:Sterile container Temperature:Sample must be maintained at body temperature	89321
LRH	SEROTON	SEROTONIN	
LRH	SEXHORM	SEX HORMONE BINDING GLOBULIN	
LRH	SIROLIM	SIROLIMUS	
LRH	SMITH	SM (SMITH) IGG AUTOANTIBODIES	
LRH	SMEAR	SMEAR REVIEW Container:Lavender-top (EDTA) tube Temperature:Ambient	85060