



Online Test Menu

Laboratory Catalog

LAB	CODE	TEST NAME	CPT4 CODE
Mission Health	LACOLC	Lacosamide Level LC Place Of Service:LabCorp	
Mission Health	LDH	Lactate Dehydrogenase via Ortho Microslide, 1LAB.CHM.5027.00 Container:4 mL SST or 4 mL PST Temperature:RT	83615
Mission Health	LA	Lactic Acid Via Ortho Microslide 1LAB.CHM.5026.00 Container:4 mL Gray Sodium Flouride Tube Temperature:Heparinized plasma on ice	83605
Mission Health	LAMLC	Lamotrigine (Lamictal) Level, Serum LC Place Of Service:LabCorp	
Mission Health	BLD	LDH, Body Fluid Container:Collect 1 mL of body fluid in STERILE container: Pleural Fluid Peritoneal Fluid Synovial Fluid Thoracentesis Fluid	83615
Mission Health	CLD	LDH, CSF Container:1.5 mL CSF	83615
Mission Health	LDL DLC	LDL Cholesterol (Direct) LC Place Of Service:LabCorp	
Mission Health	LEADLC	Lead, Blood (Adult) LC Place Of Service:LabCorp	
Mission Health	LEADPLC	Lead, Blood (Peds) Capillary LC Place Of Service:LabCorp	
Mission Health	LEADP VLC	Lead, Blood (Peds) Venous LC Place Of Service:LabCorp	
Mission Health	ULGAG	Legionella Antigen, Urine Container:5 mL Urine in STERILE CUP	87449
Mission Health	LPNELC	Legionella pneumophila DFA LC Place Of Service:LabCorp	
Mission Health	LEPSLC	Leptin, Serum LC Place Of Service:LabCorp	
Mission Health	PATHFLOW	Leukemia/Lymphoma Panel Container:Any of the following are acceptable specimen types: 3 mL Whole Blood in a 4 mL Green Sodium Heparin Tube 3 mL Body Fluid Sterile Specimen Container 3 mL Bone Marrow in a 4 mL Green Sodium Heparin Tube Other: Sterile Specimen Container	88184 88185 each additional marker; 88187 or 88188 or 88189
Mission Health	LYSEN	Leukocyte Lysosomal Enzyme Container:5 mL Green Sodium Heparin Tube Temperature:Critical Ambient Place Of Service:Jefferson Neurogenetics	82657