

# PMG Laboratory Spokane

## Test Directory

| LAB    | CODE     | TEST NAME                                                                                                                                                              | CPT4 CODE   |
|--------|----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------|
| PMGLAB | LAB25514 | EDOXABAN LEVEL<br>Place Of Service:LabCorp                                                                                                                             |             |
| PMGLAB | LAB25368 | EGFR ONCOLOGY, FISH<br>Place Of Service:LabCorp                                                                                                                        |             |
| PMGLAB | LAB16    | ELECTROLYTE PANEL<br>Container:Red top tube; SST gel barrier tube; PST gel barrier tube (green lithium heparin)<br>Temperature:Refrigerate<br>Place Of Service:PMG-PMP | 80051       |
| PMGLAB | LAB774   | ENDOMYSIAL AB, IGA<br>Place Of Service:LabCorp                                                                                                                         |             |
| PMGLAB | LAB1333  | ENTEROVIRUS, RNA, NAAT<br>Place Of Service:LabCorp                                                                                                                     |             |
| PMGLAB | LAB298   | EOSINOPHIL COUNT<br>Container:EDTA lavender top tube<br>Temperature:Refrigerate<br>Place Of Service:PMG-PMP                                                            | 85048       |
| PMGLAB | LAB1812  | EPSTEIN-BARR VIRUS AB, IGG AND IGM<br>Container:SST, red top tube<br>Temperature:Refrigerate<br>Place Of Service:PMG-PCS                                               | 86664 86665 |
| PMGLAB | LAB1730  | EPSTEIN-BARR VIRUS CAPSID AB, IGG<br>Place Of Service:LABCORP SEATTLE                                                                                                  |             |
| PMGLAB | LAB1731  | EPSTEIN-BARR VIRUS CAPSID AB, IGM<br>Place Of Service:LABCORP SEATTLE                                                                                                  |             |
| PMGLAB | LAB654   | EPSTEIN-BARR VIRUS EARLY ANTIGEN AB, IGG<br>Place Of Service:LABCORP SEATTLE                                                                                           |             |
| PMGLAB | LAB793   | EPSTEIN-BARR VIRUS NUCLEAR ANTIGEN AB, IGG<br>Place Of Service:LABCORP SEATTLE                                                                                         |             |
| PMGLAB | LAB863   | EPSTEIN-BARR VIRUS PANEL<br>Place Of Service:LABCORP SEATTLE                                                                                                           |             |
| PMGLAB | LAB2372  | EPSTEIN-BARR VIRUS, NAAT, QUAL<br>Place Of Service:LABCORP BURLINGTON                                                                                                  |             |
| PMGLAB | LAB1373  | EPSTEIN-BARR VIRUS, NAAT, QUANT<br>Place Of Service:LABCORP BURLINGTON                                                                                                 |             |
| PMGLAB | LAB24978 | EPSTEIN-BARR VIRUS, PLASMA, QUANT PCR<br>Place Of Service:LabCorp                                                                                                      |             |