

PMG Laboratory Spokane

Test Directory

LAB	CODE	TEST NAME	CPT4 CODE
PMGLAB	LAB347	UA DIPSTICK ONLY Container:Leakproof plastic urine container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	81003
PMGLAB	LAB12834	UA DIPSTICK WITH MICROSCOPIC IF INDICATED CULTURE IF INDICATED Container:Leakproof plastic urine container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	81003
PMGLAB	LAB1839	UA PROTEIN, DIPSTICK Container:Leakproof plastic container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	81002
PMGLAB	LAB2501	UAM REFLEX CULTURE Container:Leakproof plastic urine container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	
PMGLAB	LAB748	UREA NITROGEN, URINE Place Of Service:LabCorp	
PMGLAB	LAB23927	UREAPLASMA/MYCOPLASMA HOMINIS Place Of Service:LabCorp	
PMGLAB	LAB141	URIC ACID Container:PST gel barrier tube; SST gel barrier tube; Red top tube Temperature:Refrigerate Place Of Service:PMG-PMP	84550
PMGLAB	LAB2479	URINALYSIS WITH MICROSCOPIC IF INDICATED Container:Leakproof plastic urine container or urine transport tube Temperature:Refrigerate Place Of Service:PMG-PMP	
LABCORP	511200	UGT1A1 Irinotecan Toxicity Container:Lavender-top (EDTA) tube, yellow-top (ACD) tube OR LabCorp buccal swab kit Temperature:Maintain whole blood specimen at room temperature for 7 days or refrigerated for 1 month, frozen for 2 years. Maintain buccal swabs at room temperature for 6 months. **Do NOT Freeze glass tubes**	81350
LABCORP	806050	Ulocladium chartartum IgE Container:Gel-barrier tube Temperature:Maintain specimen at room temperature	86003
LABCORP	470115	Uniparental Disomy (UPD) Parent, DNA Analysis Container:Lavender-Top (EDTA) Tube, Yellow-top (ACD) tube Temperature:Room temperature. Maintain specimen at room temperature. Do not freeze.	NOCPT
LABCORP	470074	Uniparental Disomy (UPD) Proband, DNA Analysis Container:Lavender-Top (EDTA) Tube, yellow-top (ACD) tube; sterile plastic conical tube or two confluent T-25 flasks for fetal testing Temperature:Room temperature. Maintain specimen at room temperature. Do not freeze.	NOCPT
LABCORP	813253	UPG III Synthase,Erythrocytes Container:Green-top (heparin) tube Temperature:Refrigerate	82657
LABCORP	008342	Upper Respiratory Culture, Routine Container:Bacterial culture transport swab Temperature:Maintain specimen at room temperature.	87070
LABCORP	808567	Uranium, Urine Container:Metal-free plastic urine container or metal-free transfer tube Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	83018