

Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB24739	PAIN MANAGEMENT BENZODIAZEPINE REFLEX Place Of Service:LabCorp	
SPH	LAB24791	PAIN MANAGEMENT PANEL 3 REFLEX Place Of Service:LabCorp	
SPH	LAB24783	PAIN MANAGEMENT PANEL 4 REFLEX Place Of Service:LabCorp	
SPH	LAB979	PANCREATIC ELASTASE, FECAL Place Of Service:LabCorp	
SPH	LAB16084	PANCREATIC POLYPEPTIDE Place Of Service:LabCorp	
SPH	LAB1790	Pap Smear Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141
SPH	LAB12010	PAP, LIQUID BASED Place Of Service:LABCORP SEATTLE WA CYTO	
SPH	LAB7088	Pap, Liquid Based, and CT/NG Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141 87491 87591
SPH	LAB7098	Pap, Liquid Based, and HPV, High Risk Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141 87624
SPH	LAB13425	PAP, LIQUID BASED, AND HPV, HIGH RISK Place Of Service:LABCORP SEATTLE WA CYTO	
SPH	LAB24969	Pap, Liquid Based, and HPV, Reflex, Genotype 16,18,45 Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141 87625
SPH	LAB12009	Pap, Liquid Based, Reflex HPV Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141 87624
SPH	LAB7087	Pap, Liquid Based, Reflex HPV if ASCUS, and CT/NG Container:ThinPrep® collection vial, broom-type collection device with detachable head, Cervix brush with snap off head. Temperature:Room Temperature Place Of Service:Incyte Diagnostics (Spokane)	88175 88142 88141 87624 87491 87591
SPH	LAB23708	PAP, LIQUID BASED, REFLEX HPV, REFLEX GENOTYPE 16,18 Place Of Service:LABCORP SEATTLE WA CYTO	
SPH	LAB24785	PARACOCCIDIODES BRASILIENSIS TOTAL ANTIBODY, IGG AND IGM Place Of Service:LabCorp	